

COLPOSCOPE HD

SERIES C-100

- Outstanding image quality • Perfect HD 1280X720p color definition • Optical zoom
- LED high power illumination • Great maneuverability



COLPOSCOPE HD

SERIES C-100



Optional:

- Monitor stand pole
- Endodigi software capturing system.
- JPG output



TECHNICAL SPECIFICATIONS

SERIES	LED LIGHT SOURCE
OUTPUTS	2x HDMI
MAGNIFICATION	progressive magnification. Optical zoom up to 50x and digital zoom up to 128x.
FOCUS	Automatic and manual
RESOLUTION	1280x720P
DEPTH OF FIELD	250-300
FIELD OF VIEW	8-140 mm
DISTANCE	200-320mm.
FILTER	optical green
LIGHT SOURCE	Led 6500 K
DURABILITY	50000 hs.
VOLTAGE	100-240 volts.
FREQUENCY	50/60hz
WORKING TEMPERATURE	5-40 degree Celsius
HUMIDITY	less than 80% without condensation
ATMOSPHERIC PRESSURE	70-106k

- OSD (on screen display) with indication of real magnification.
- Zoom and focus operated by 4 commands joystick.
- Still function.
- WB (white balance) push bottom control.
- Perfect color discrimination due to CMOS HD technology.
- Floor stand system with 4 antistatic wheels with break.
- Ratched up down pole with double arm for easier and steady positioning.

Rev.31.03.14



www.ecleris.com

ECLERIS
INNOVATION VALUE