


E.N.T.
RIGID AND FLEXIBLE
ENDOSCOPES

ENDOSCOPES
SINUSCOPES / LARYNGOSCOPES / OTOSCOPE



ECLERIS
INNOVATION VALUE

We offer exclusive design,
excellent quality and cost
effective endoscopes.
Ecleris innovation value

ENDOSCOPES
SINUSCOPES / LARYNGOSCOPES / OTOSCOPE

Outstanding image quality.
Exceptional luminosity and brightness.
Storz / Wolf and ACMI connections.



ENDOSCOPES	
DIAMETER	4.5mm
WORKING LENGTH	Approx. 140mm
LOCKING DEVICE	Storz Type
DIRECTION OF VIEW	0° 30° 70°



ENDOSCOPES	
DIAMETER	4.5mm
WORKING LENGTH	Approx. 170mm
LOCKING DEVICE	Storz Type
DIRECTION OF VIEW	0° 30° 40° 70°



OTOSCOPE	
DIAMETER	3.5mm
WORKING LENGTH	Approx. 100mm
LOCKING DEVICE	Storz Type
DIRECTION OF VIEW	0°



LARYNGOSCOPE	
DIAMETER	5.5mm
WORKING LENGTH	Approx. 170mm
LOCKING DEVICE	-
DIRECTION OF VIEW	70° 70°



FLEXIBLE SCOPE	
FIELD OF VIEW	80°
DEPTH OF FIELD	2.5 - 50 mm
OUTER DIAMETER INSERTION TUBE	3.8 mm
OUTER DIAMETER OPTIC END	2.4 mm
WORKING LENGTH	300 mm
BENDING ANGLE OF DISTAL TIP	up 130° down 130°

Ref. 001171