

PROCAM HD

E.N.T
HIGH DEFINITION CAMERA

EXCELLENT COLOR REPRODUCTION

• GREATER CLARITY

• SHARP IMAGES



• IMPROVED
DEPTH PERCEPTION

• ERGONOMIC
CAMERA HEAD

• ON SCREEN DISPLAY

- WHITE BALANCE
- WINDOWS FUNCTION (2 LEVELS)
- ANTI - MOIRÉ (2 LEVELS)
- STROBOSCOPY SET UP

PROCAM HD



TECHNICAL SPECIFICATIONS

RESOLUTION	1280 X 720 p
VIDEO SYSTEM	progressive scan 60 / 50hz
CCD SENSOR	1/3" HD
SENSITIVITY	< 3 LUX
SETTINGS	on screen display
OUTPUTS	3 X HDMI 1 X COMPOSITE (NTSC/PAL)
POWER SUPPLY	100 – 240 VAC
DIMENSIONS (W X D X H)	300 X 290 X 95 cm.
HEAD	C – mount thread
CLASS	I
TYPE	B

DISTRIBUTED BY:

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE. FL-63001 Rev.10.08.11



COST EFFECTIVE SOLUTIONS FROM ECLERIS